



Reference: PTXBC052955

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	N-WASP; NWASP; WASPB; neural Wiskott-Aldrich syndrome protein; Wiskott-Aldrich syndrome like		

For research use only.

GENE DESCRIPTION	Wiskott-Aldrich syndrome-like
GENE ID	8976
ORIGIN SPECIES	Human
ACCESSIONS	BC052955
VECTOR	pDONR223
NCBI SYMBOL	WASL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]