



Reference: PTXBC052978

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pleckstrin homology domain-containing family H member 3; PH domain-containing family H member 3; pleckstrin homology domain containing, family H (with MyTH4 domain) member 3; pleckstrin homology, MyTH4 and FERM domain containing H3		

For research use only.

GENE DESCRIPTION	pleckstrin homology domain containing, family H (with MyTH4 domain)
GENE ID	member 3 79990
ORIGIN SPECIES	Human
ACCESSIONS	BC052978
VECTOR	pDONR223
NCBI SYMBOL	PLEKHH3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr