

SH3BGRL2-SH3 domain binding glutamic acid-rich protein like 2 Gene

Reference: PTXBC052987

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SH3 domain-binding glutamic acid-rich-like protein 2; SH3 domain binding glutamic acid-rich protein like 2; fovea-associated SH3 domain-binding protein; SH3 domain binding glutamate rich protein like 2		

For research use only.

GENE DESCRIPTION	SH3 domain binding glutamic acid-rich protein like 2
GENE ID	83699
ORIGIN SPECIES	Human
ACCESSIONS	BC052987
VECTOR	pDONR223
NCBI SYMBOL	SH3BGRL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcatccgctgtttcatcgccctctctcgggcttcgtggcgataaagaagcagcaagatgtggttagattctggaagccaacaagatagagtttgaggaggtggatatcacaatgtcagaagaacagaggcaatggatgtacaaaaacgtcccccgaaaagaaccctcagggaacccccctgccacctcagatatattaatggcgaccgatactgtggagattatgacagtttttgaatccaaggaagcaacacagctctttcattttaggcctgaaaccacgggttgcatcaaggcagaaccttag