



Reference: PTXBC053325

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf115; vertnin; vertebrae development homolog (pig); vertebrae development associated		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 115
GENE ID	55237
ORIGIN SPECIES	Human
ACCESSIONS	BC053325
VECTOR	pDONR223
NCBI SYMBOL	C14orf115
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr