



Reference: PTXBC053347

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S-MAG; GMA; SIGLEC-4A; SIGLEC4A; SPG75; myelin-associated glycoprotein; sialic acid binding Ig-like lectin 4A; sialic acid-binding immunoglobulin-like lectin 4A		

For research use only.

GENE DESCRIPTION	myelin associated glycoprotein
GENE ID	4099
ORIGIN SPECIES	Human
ACCESSIONS	BC053347
VECTOR	pDONR223
NCBI SYMBOL	MAG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr