



Reference: PTXBC053349

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FRTS2; HCINF2; NAPI-3; NPHLOP1; NPT2; NPTIIa; SLC11; SLC17A2; sodium-dependent phosphate transport protein 2A; Na(+)-dependent phosphate cotransporter 2A; Na(+)/Pi cotransporter 2A; Na+-phosphate cotransporter type II; naPi-2a; renal sodium-dependent phos		

For research use only.

GENE DESCRIPTION	solute carrier family 34 (sodium phosphate), member 1
GENE ID	6569
ORIGIN SPECIES	Human
ACCESSIONS	BC053349
VECTOR	pDONR223
NCBI SYMBOL	SLC34A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr