

SPINLW1-serine peptidase inhibitor-like, with Kunitz and WAP domains 1 (eppin) Gene

Reference: PTXBC053369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPINLW1; CT71; CT72; WAP7; WFDC7; dJ461P17.2; WAP four-disulfide core domain protein 7; cancer/testis antigen 71; cancer/testis antigen 72; epididymal protease inhibitor; protease inhibitor WAP7; serine peptidase inhibitor-like, with Kunitz and WAP domain		

For research use only.

GENE DESCRIPTION	serine peptidase inhibitor-like, with Kunitz and WAP domains 1 (eppin)
GENE ID	57119
ORIGIN SPECIES	Human
ACCESSIONS	BC053369
VECTOR	pDONR223
NCBI SYMBOL	SPINLW1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatctctggactttgagcctcctggctatcctcttagcgaatgtccaggacctggctgactgattggtattccaggagatgtccaaaatcagagaagaatgtgaattccaagaaaggatgtgtgtacaaggacagacaatgccaggacaacaagaagtgtgtgtcttcagctcggaaaaaatgttagatcacaacaagatgtatgcgaatgccaaaagaactgcccctgctgcttatttctcattggtgtatgacaagaaagataactgtctcatgtttgtctatggtggctgccagggaacaataaacaactccaatccaagccaactgcctgaacacctcaagaataaacgcttccctga