

## TRIM16-tripartite motif-containing 16 Gene

Reference: PTXBC053514

|                |                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | EBBP; tripartite motif-containing protein 16; estrogen-responsive B box protein; tripartite motif containing 16 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | tripartite motif-containing 16                                                                   |
| <b>GENE ID</b>          | 10626                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC053514                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | TRIM16                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggctgagttggatcattggtccaggccactgccagcccccagccctctcagccagactctgggtcaccagccagattctgggtcagccagccagtggaagaagaggacgtgggtcctcggagaagctggcaggagacggaggaacagga  
cagcgactctgcagagcaggggatctgctggtgaggggaaagaggtctgtgtgactctgctgtgatgacaccagaagagtgaggcagtgaaagtcctgtctaacctgcagtggaattactgtgaaagacactgcagccgcacatcagtggaacatcaaacgcaagccactgct  
gaccgagccagtggaaggaccacaactggcgatactgacctgccaccagccactgtctgctctgctgacctgatcagcagtgcatctgccaggactgttgccaggagcacagtgccacacatagctctccctggatgcagcccgagggacagggaggtgaactccagtg  
accagtttagacttggagcggaaactcaagttgaatgaaatgccatctccaggctccaggctcaacaaaagtctgtctggtgtcgggtgtcagaggtcacaagcgggtggctgaaatgcagtttggggaactccttgcctgtgtgaggaagcccgccaatgtgatgctctcttagagga  
gaaggagcaagctgcctgagccagcccaagcgtatcaagggccacctggagtacaggagtgccgagatggagaagagcaggagctggagagggatggcggccatcagcaacactgtccagttcttgaggagtgactgcaagtttaagaacactgaagacatcaccttcccta  
gtgtttacgtaggctgaaggataaactctcggcatccgaaagtatcacggaaatccactgtacacttaaccagttgtctggagaactataagaaaaagctccaggagtttccaaggaagaggagtgatgacatcagaactcaagtgctgctggtgttcagcgcgaatattggacttccaa  
cctgagcccgagcaccaggggaacagttcctccaatgctgtagatcacgtttgacccggacacagcacaaagtatctccggctgcaggaggagaaacccaaggtcaccacacacacgcccctgggagcatccctaccggacgtcccccagcaggttctgctgctgctgctgctg  
ctgtccagcagagctgtacctgcacaggtactattttgaggtggagatcttccggcagcagcactatgttggcctgacctgcaaaagcagtcgacggaaagggaggagcgcacaggttgcatttccggaacaaacttctctggagcctccaatgggaacgggaaggaggtcacggcc  
tggtacagtgacatggagacccactcaagctggcccttccggaggtcggggtctatatcgacttccgggagggatccttctctatggcgtagatgataccatgactctgttcacaaagtttgcctgcaaatitcagaacacagtgatgctgcttctggttccaagaaggaaa  
acgcatccggttagatctgggagaggaaaccgagaagccagcaccgtccttgggtgactgctccctag