

ATP9B-ATPase, class II, type 9B Gene

Reference: PTXBC053561

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATPASEP; ATP1IB; HUSSY-20; NEO1L; hMMR1; ATPase type IV, phospholipid transporting (P-type); ATPase, class II, type 9B; macrophage MHC receptor 1; ATPase phospholipid transporting 9B (putative)		

For research use only.

GENE DESCRIPTION	ATPase, class II, type 9B
GENE ID	374868
ORIGIN SPECIES	Human
ACCESSIONS	BC053561
VECTOR	pDONR223
NCBI SYMBOL	ATP9B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccactaatgatgtctgaagaagccttgagaatgaggaaagtgattaccacacctaccacgagccaggataatgcaaggaaaaggaggactggagtggtttgtctgtgtagggctggaagttcctctgtaccagttgctgtggttgctgataaatattgtcgaagaagaagagctga
aagctgcacagtatggcttgatgtctgaaaagtgtgagaaaaacatccaggaaattctataaaaaatcaaaaatacaatgtgtttacctttatactggggtttgtatgaacaattcaagttttctgaatctctattttctagtaaatatcctgctcacagttgtaccagcattgaaaataggcta
tctctacacctactgggctcctctgatgacatga