

NXPH4-neurexophilin 4 Gene

Reference: PTXBC053581

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	alternative protein NXPH4; NPH4; neurexophilin-4; neurexophilin 4		
For research use only.			

GENE DESCRIPTION	neurexophilin 4
GENE ID	11247
ORIGIN SPECIES	Human
ACCESSIONS	BC053581
VECTOR	pDONR223
NCBI SYMBOL	NXPH4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcggctgtcccggaatggtctcttctgtcttggccgtggctcttaggaaggccgtcagtgccagataccagatccggaaggccgagtcctgggctgcgccccgcgcggcgagcgggtgccccggccagcagctccagagccaaggcttcggacggccta
ggcgtggggccgcgcctggagctgggctggccgaccaaccacacggggcgctggccgggagggcgagccggggcgttggccgcgcagcgcaccaagaggagccgtccatcaaggcgccgcgcgcaaaaagatctcggctgggggacttctacttctgggtgcatac
cctcaagtttctgtgtggtgaccggcaagatcgtggaccatgtgaacgttaccttcagtggtgtatctccgcaaacctgtccagcctgggcaacctcagtgtagcatcgtgcccctccaagcgtgtcagttcggaggagctctggctgccccggcctgtccccaccctctgcagtc
tacgtcgcgccctggaggggggtcttctggtgggcccccgctggggatggcagcagcagcggcgccggggcgttggggcgtccctcggggcgcaactggcgggcgccgttggggcgcggttggggagtgccctggggcgaagagtcacgcgcttcaattgccacgtgga
gtatgagaagacaacccgcgcgcaagcaccgacgtgctgtacgatccatcgaggtgtgtcttcaccgagcacacgcagagccagggccgtgctgtgccaagccctcaaaagtcgtgtatcttctcgtcttctcagctttgactacaaactggtgcagaaggtgtgccag
actataacttccagatgagcaccctacttccgatag