



Reference: PTXBC053600

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane and coiled-coil domain-containing protein 4; transmembrane and coiled-coil domains 4		

For research use only.

GENE DESCRIPTION	transmembrane and coiled-coil domains 4
GENE ID	255104
ORIGIN SPECIES	Human
ACCESSIONS	BC053600
VECTOR	pDONR223
NCBI SYMBOL	TMCO4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]