

## C15orf38-chromosome 15 open reading frame 38 Gene

Reference: PTXBC053602

|                |                                                                                                            |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                        | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                       | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | UPF0552 protein C15orf38; C15orf38; arp2/3 inhibition protein; actin-related protein 2/3 complex inhibitor |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 15 open reading frame 38                                                               |
| <b>GENE ID</b>          | 348110                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC053602                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C15orf38                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgagccgcatctaccagacgagcgctccggaacaggcgggtgcagagcgtccggctccaggggctgggacccgcccaccaggggggaaatggtgctctgctggaggagaaatgatgatgtatctcggcacagcatcttgacactcatggcagggaaggagcgt  
actacgtgctgtatctcggccagtcacatccatcgcgttaattcgacgccaagggaatgaaatcgaagcccaactcagcgcaccaggaagggtgaacacgggcttctcatgtctctacaagggtggaagccaagggggacactgacaggtcacgcccaggcgtgaagg  
ggctggtcaacaagccagagctgctcgcgtgacagagagcctcaccgccaccacagctggtgctggtgagtcgcccagtcagagatggaggtgatgaactcagctggggccgggggtacggctgaagactcggggcgtggtccctcctggatcattggccaactgagg  
ctggaacagtgaccaagtgttaattcactggtgatgaaagacaggggcatcctggacagacaacatcatggccaaaagtgttcgaaggggctgcagcggagatccgagagcaggggggatgggagaggacagagagtggtgatgactga