

PRKAB2-protein kinase, AMP-activated, beta 2 non-catalytic subunit Gene

Reference: PTXBC053610

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	5'-AMP-activated protein kinase subunit beta-2; 5'-AMP-activated protein kinase, beta-2 subunit; AMP-activated protein kinase beta 2 non-catalytic subunit; AMPK beta 2; AMPK beta-2 chain; AMPK subunit beta-2; protein kinase, AMP-activated, beta 2 non-cata		

For research use only.

GENE DESCRIPTION	protein kinase, AMP-activated, beta 2 non-catalytic subunit
GENE ID	5565
ORIGIN SPECIES	Human
ACCESSIONS	BC053610
VECTOR	pDONR223
NCBI SYMBOL	PRKAB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaacaccaccagcagccgggtgtccggggagcgcacagggcgtgcacgctccgagggcgacggcggccatgccccggggaaggagcacaagatcatgggtgggagtagcagcagccccagcgtgttcagcctccctgactccaagctccctggggacaaag
agtttgtatcatggcagcaggtttggaggactccgtaaaagccacacagcagggccggccactgttatccgctgggtctgaaggaggcaaggaggtcttcatctctgggtcttcaacaattggagcaccagattccactgattaagagccataatgactttgtccatcctggacctccc
tgaggagagcaccatacaagttcttggatggacagtggttcatgatccatcagagcctgtgttaccagtcagcttggcacaattacaattgatccatgtcaagaatctgattttgaggtgttcgatgctttaaagtttagattctatgaaagtctgagacatcttagagacattcc
agctcaccggcggccttatggtaagaaatgtatcggttcgatctgaggaagattcaaatcccacccatcctctctcatctcaagttattcttaacaaagacactaatatttcttgaccagccttactccctgagcccaaccatgttatgtgaacctctctatgcattgtccat
taaggacagtgatgttccttagcgcgaacccatcgctacaagaagaagtattgtactactctgctatacaagccatttga