

NRARP-NOTCH-regulated ankyrin repeat protein Gene

Reference: PTXBC053618

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	notch-regulated ankyrin repeat-containing protein; NOTCH-regulated ankyrin repeat protein		

For research use only.

GENE DESCRIPTION	NOTCH-regulated ankyrin repeat protein
GENE ID	441478
ORIGIN SPECIES	Human
ACCESSIONS	BC053618
VECTOR	pDONR223
NCBI SYMBOL	NRARP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccaggccgagctgtccacctgctccgcccagacgcagcgcattccaggagctgtgcgaagggaacacgcaggagctgcagtcgctgctgcagaacatgaccaactgcgaggtcaacgtgaactcgttgggcccagggccagacggcgctgcaccagtcgg
tcacgcagggcaacctggagctcgtgaagctgctgtaagttcggccgacatccgctggccaaccgcgacgctggagcgcgctgcacatgccgcgttcggtggccaccaggacatcgtgctctatctcatccaaggcgaagtacggccagcgccggtga