

C5orf24-chromosome 5 open reading frame 24 Gene

Reference: PTXBC053677

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0461 protein C5orf24; chromosome 5 open reading frame 24		
For research use only.			

GENE DESCRIPTION	chromosome 5 open reading frame 24
GENE ID	134553
ORIGIN SPECIES	Human
ACCESSIONS	BC053677
VECTOR	pDONR223
NCBI SYMBOL	C5orf24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgcctcgtgtgccagcagtaatccagcttctgtggcctggcaagccttcctgcctcaatgaagatgccatgagagctgtgatcagttgacatatattcctccagcaagcaaatagccacacagtcacccaacaaatggtttgtcagaggcaagaccattaaatgaaa
cacactgcagactacaagtggcagaagcatagaaataaaagatgaactaaagaaaaagaagaatctcaaccgatctgtaagcgtggccggccttcgggaaccaccaaatcagcaggataccggaccagcacaggcagaccctgggaaccaccaagcagctggatttaagacaa
gtccaggcagacctttggggacaactaaagctgctgggatacaaaagtcagcccagcagacctccaggtagcattaaagctctatccgctcttgcgatcttgggtatggctgtggcactgctgttttccttacctatgatgatggcagagcagttcatgggtagaggaaactagcagtg
aagtcacacccaatgagtg