

KLRC1-killer cell lectin-like receptor subfamily C, member 1 Gene

Reference: PTXBC053840

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD159A; NKG2; NKG2-A/NKG2-B type II integral membrane protein; C-lectin type II protein; CD159 antigen-like family member A; NK cell receptor A; NKG2-1/B activating NK receptor; NKG2-A/B type II integral membrane protein; NKG2-A/B-activating NK receptor;		

For research use only.

GENE DESCRIPTION	killer cell lectin-like receptor subfamily C, member 1
GENE ID	3821
ORIGIN SPECIES	Human
ACCESSIONS	BC053840
VECTOR	pDONR223
NCBI SYMBOL	KLRC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggataaccaaggagtaactctactcagacctgaatctgcccccaacccaaggaggcagcaacgaaacctaaggcaataaaagctccattttagcaactgaacaggaaataacctatgcggaattaaacctcaaaaagctctcaggatttcaagggaatgacaaaacctatcactg
caagatttaccatcagctccagagaagctcattgttgggatctgggaattatctgtcttcttaattggcctctgtggaacgatagttgttattccctctacattaatacagaggcacaacaattctccctgaatacaagaactcagaagcagctcattgtggccattgtcctgaggagtgatt
acatattccaacagttgttactacattgtaaggaaagaagaacttgggaagagatttgcctggcctgtacttcgaagaactccagctcgttctctatagataatgaagaagaatgaattctgtccatcattcaccatcctcatggattggtgtgttcgtaacagcagtcacatccatgggtg
acaatgaatggttgccttcaacatgagataaaagactcagataatgctgaacttaactgtgcagtgtctacaagtaaatcgaactaaatcagcccagtggtgatctcaataatatacattgtaagcataagccttag