

APOBEC3B-apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3B Gene

Reference: PTXBC053859

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A3B; APOBEC1L; ARCD3; ARP4; DJ742C19.2; PHRBNL; bK150C2.2; DNA dC->dU-editing enzyme APOBEC-3B; apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3B; cytidine deaminase; phorbol 2; phorbol 3; phorbol-1-related protein; phorbol-2/3;		

For research use only.

GENE DESCRIPTION	apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3B
GENE ID	9582
ORIGIN SPECIES	Human
ACCESSIONS	BC053859
VECTOR	pDONR223
NCBI SYMBOL	APOBEC3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatccacagatcagaaatccgatggagcggatgacgagacacattctacgacaacttgaaaacgaacccatcctctatggtcggagctacacttgctgtgctatgaagtgaataaagagggccgctcaaatctcttgggacacaggggtcttcgagccaggtgtatttc
gagcctcagtagccagcagaaatgcttctctcttgggtctgtggcaaccagctgctgttccagatcacctggtttgtatcctggacccccctggcggactgtgtggcgaagctggcgaattctgtctgagcaccccaatgtcaccctgacctctccgccccgct
ctactactctgggaagagattaccgaagggcgctctgagcgtgagtcaggcaggagcccgctgaagatcatggactatgaagaatttgcatactgtgggaaaacttgggtacaatgaagtcagcaattcatgcttggtaaaattcgatgaaattatgattctgcaccgcac
gctaaaggagatttcaggctaagaatctctctgtggttttacggcccatgcggagctgcctcttggacgtggttctcttgcagttggacccggccagatctacaggtgcttgggtcatctctgagccctgcttctctggggtgtgccccgggaagtgcgtggttcttc
aggagaacacacacgtga