

CPNE5-copine V Gene

Reference: PTXBC053872

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COPN5; CPN5; copine-5; copine V; copine 5		

For research use only.

GENE DESCRIPTION	copine V
GENE ID	57699
ORIGIN SPECIES	Human
ACCESSIONS	BC053872
VECTOR	pDONR223
NCBI SYMBOL	CPNE5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaaaaaatacgtgaattctggcacagtcaccctgtttcctttgtgtggagtcagagtgacacctcttgactacatcaaaggaggaccagatcaacttcactgtggccattgatttcactgcctcaatgggaacccctcacagtcacatccctgcactacatgagcccta
ccagctgaacgcctacgcgctggcgtgactgcccgtggagagatcatccagcactacgacagtgacaagatgttccctgcccgtggccttccgggccaagctgcccccgatggcagagtgctccacgagttcccaatgaatggcaaccaggagaacccctcatgctgcgcatcgac
ggcatcctggaggcctaccaccgcagcctgacgactgtgacgtgtacggcccccaactttgccccgtgtgacccacgtggccaggaatgcagcggcgtgaggaatggctcccaatgactggtgctgctcatcattactgatgggcatctggacatggcgagaccaagga
ggccattgtcaacgtgccaagctcccatgtccatcattatctgctggcggccagagagttcgacgcatggtggagctggatggcgacgacgtgcgcatctctccgggggaagctggtgaacgcgacatcgtccagtgtgtacccttcgggactacgtggaccgcacag
gcaaccacgtgctgagcatggccgcctggcccgagacgtgctggcagagatccctgaccaactggtgtcctacatgaaggcacagggcattgcggcgtgccccaccgcagcaccacccactgcccctgcagtcgccagcccgaogccccctgctccccctgcacacg
cacatctga