

CCDC117-coiled-coil domain containing 117 Gene

Reference: PTXBC053874

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dJ366L4.1; coiled-coil domain-containing protein 117; coiled-coil domain containing 117		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 117
GENE ID	150275
ORIGIN SPECIES	Human
ACCESSIONS	BC053874
VECTOR	pDONR223
NCBI SYMBOL	CCDC117
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcgctcggccgctcagcgccctcctctgagcggcgctcggacttcctcagcggccgagcggcctcctccggcgggcctcggcggggctgacggcgccgagttggcccgcgccgggacctcgcgcagtccttagcagtcgggagtgaggcgcgcgacgtgttctgttctgactgtaaaaagaacacacagcagaggaggagggatgatgtgtccagttaagaaagaaagataactgaagcagagctctgtgctggtcctaatgactggattctttgtgcacatcaggatgtagagggcatggagtaaatccagtgtagtgcccttcatactggatattagatgttattgtgaagaaatggatcagacaactggagaaccacagtggtgaagttgccgaaggagcttcaggagattgaggacaggataattgatgaagatgaagaagttgaagctgacagaaatgttaacctctcccagcttctccttctgataccatgaaaacaggtttgaagagggaattgatgaagttttacaagaataatgattgagctctatgagccgtcctccatggagcttctctctggaacccctccctgaactccttctgataagccaagccatccttaataactaagaactatacaggagagagccaagctaagcatgtagctgctgcactgcctccctcagagaactgaactgttttcggaacctcgccaacaggagatgtctcttataatagtttggagacagctactagcacagaagaagatggaacttag