



Reference: PTXBC053877

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTL3; choline transporter-like protein 3; solute carrier family 44 member 3		
For research use only.			

GENE DESCRIPTION	solute carrier family 44, member 3
GENE ID	126969
ORIGIN SPECIES	Human
ACCESSIONS	BC053877
VECTOR	pDONR223
NCBI SYMBOL	SLC44A3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]