



Reference: PTXBC053881

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD316; CD81P3; EWI-2; EWI2; KCT-4; LIR-D1; PGRL; immunoglobulin superfamily member 8; CD81 partner 3; glu-Trp-Ile EWI motif-containing protein 2; keratinocytes-associated transmembrane protein 4; prostaglandin regulatory-like protein		

For research use only.

GENE DESCRIPTION	immunoglobulin superfamily, member 8
GENE ID	93185
ORIGIN SPECIES	Human
ACCESSIONS	BC053881
VECTOR	pDONR223
NCBI SYMBOL	IGSF8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]