

HOXC9-homeobox C9 Gene

Reference: PTXBC053894

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOX3; HOX3B; homeobox protein Hox-C9; homeo box C9; homeobox protein Hox-3B; homeobox C9		

For research use only.

GENE DESCRIPTION	homeobox C9
GENE ID	3225
ORIGIN SPECIES	Human
ACCESSIONS	BC053894
VECTOR	pDONR223
NCBI SYMBOL	HOXC9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggcgacggggccatcagtaactattacgttgactcgtcatctctcacgacaatgaagacctctagcgtccaggtttccggccaccggggctcatcccgccgcccagaccagcgggttggtgccgactgtagcgaatttcctctgtagcttcgcgccaagccggcag
tgttcagcacgtcgtggcgcccggtgccctctcagtcgtccgtgtatatacccgtaaggccccagccccacctcggcgccgacacgcgtacatcggaacttggtcgaagccgtgtccggcgccgtctctctcccgagctccggccggggcgctcactacgcccctaaagccg
gacgctacccccggcgccgctggcgtgactgcggccagggggggcgagctacccggactacatgtacggctcggcggggagctgcgcgacgcgccccgcagacactgccccggcgagcgcgtccggcgagcaagcacaagaggagaaggctgacct
ggacccagcaacccgtggccaactggattcacgcccgtccaggaagaagcgtgccctacaccaagtaccagacgtggaactggagaaggagttctctcaatatgtatitaaaccaggaccgtcgtatgaggtggcccggttctcaatcaccgagcggcaggtca
aaatctggttcagaatcgaaggatgaagaatgaaaaatgaataaagagaaaaccgacaaggagcagtcctaa