

IRS1-insulin receptor substrate 1 Gene

Reference: PTXBC053895

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HIRS-1; insulin receptor substrate 1		
For research use only.			

GENE DESCRIPTION	insulin receptor substrate 1
GENE ID	3667
ORIGIN SPECIES	Human
ACCESSIONS	BC053895
VECTOR	pDONR223
NCBI SYMBOL	IRS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr