

MYOG-myogenin (myogenic factor 4) Gene

Reference: PTXBC053899

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MYF4; bHLHc3; myf-4; class C basic helix-loop-helix protein 3; myogenic factor 4		
For research use only.			

GENE DESCRIPTION	myogenin (myogenic factor 4)
GENE ID	4656
ORIGIN SPECIES	Human
ACCESSIONS	BC053899
VECTOR	pDONR223
NCBI SYMBOL	MYOG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctgtatgagacatccccctacttctaccaggaaccccgcttctatgatgggaaactacctgcctgtccacctcagggttcgaaccacaggtctacgagcggagctcaccctgagccccgagccccagggcccttgaggacaagggctggggacccccgag
cactgtccagggcagtgccctgccgtggggcgtgtaaggtgtgtaagaggaagtggtgtccgtggaccggcgccggcgccacactgagggagaagcgcaggtctcaagaagtgatgagccttcgagggcctgaagagaagcaccctgtcaacccaaccagcggctgcc
aaggtggagatcctgcgcagtgccatccatgacatcagcgcctccagggccctgcctcagctccctcaaccaggaggagcgtgacctccgtaccggggcggggcgcccccagccaggggtgccagcgaatgcagctctcacagcgcctcctgcagtcagagtgaggcagtg
cactggagttcagcgccaacccaggggatcctgtcagcgtgacccctacagatgccacaacctgcactccctcacctccatcgtggacagcatcacatgggaagatgtgtctgtggcctccagatgaacatgcccaactga