



Reference: PTXBC053989

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BFP; ZNF179; RING finger protein 112; brain finger protein; neurolastin; zinc finger protein 179		

For research use only.

GENE DESCRIPTION	ring finger protein 112
GENE ID	7732
ORIGIN SPECIES	Human
ACCESSIONS	BC053989
VECTOR	pDONR223
NCBI SYMBOL	RNF112
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]