

FRG1-FSHD region gene 1 Gene

Reference: PTXBC053997

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FRG1; FRG1A; FSG1; FSHD region gene 1 protein; facioscapulohumeral muscular dystrophy region gene-1; FSHD region gene 1		

For research use only.

GENE DESCRIPTION	FSHD region gene 1
GENE ID	2483
ORIGIN SPECIES	Human
ACCESSIONS	BC053997
VECTOR	pDONR223
NCBI SYMBOL	FRG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgagtagctactctgaagctacccaagctcgtctcaagggaaccaagacgaagagtaagaaagaaagagcaagataagaaagaaagagagaagaatgaagaaaccagcttgatattgttgaatctggtggacagtaacaaacttggtaaatcaggaaccata
gccattgaaatggataagggaaacctatatacatgcactcgacaatggtctttttaccctgggagctccacaaagaagttgatgaggccctagtcctccagagcagttttacggctgtcaaatatctgattccagaatcgcctgaaatctggctatggaaaatcttggataaattcagatg
gactgtgttggcgctcagatgcaattggaccaagagaacaatgggaaccagctcttcaaaatgggaaatggccttggcctcaaatagctgctttattagatgcaatgaagcagggacatagaagcaaaagtaaaacagcaggagaagaagaatgatcaagattagatcctgtg
ctgaaagagaaacaaagaaagatgacattccagaagaagacaaaggaatgtaaaacaatgtgaatcaattatgaaagaaattcagagctccagaccacaaactaaataagtaaaagagacagtaaaattctaaaaaggctcggaaagatggattttgcatgagacgcttc
tgacagggagagccaaattgaaagccgacagatactgcaagtga