

HIST1H4C-histone cluster 1, H4c Gene

Reference: PTXBC054014

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H4/g; H4FG; dJ221C16.1; histone H4; H4 histone family, member G; histone 1, H4c; histone cluster 1, H4c; histone cluster 1 H4 family member c		

For research use only.

GENE DESCRIPTION	histone cluster 1, H4c
GENE ID	8364
ORIGIN SPECIES	Human
ACCESSIONS	BC054014
VECTOR	pDONR223
NCBI SYMBOL	HIST1H4C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgtcgcgcgcaaggcggaaaaaggcttgggaagggtgtgtgtaagcgccatcgtaaggtgctcgggataacatccaggcattacaaaacggctattcgcggtttggctcggcgcggtggcgtcaagcgcatctccggtcttatctatgaggagactcgaggtgtgcttaaggtttcttagagaacgttatcagagacgccgtcacctatacgggacacgcaagcgcaaaactgtcacagccatggatgtatgtatgcctaaaacgtcagggcgcaactctgtatggcttcggcggtga