

PCBD2-pterin-4 alpha-carbinolamine dehydratase/dimerization cofactor of hepatocyte nuclear factor 1 alpha (TCF1) 2 Gene

Reference: PTXBC054021

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DCOH2; DCOHM; PHS2; pterin-4-alpha-carbinolamine dehydratase 2; 4-alpha-hydroxy-tetrahydropterin dehydratase 2; 6-pyruvoyl-tetrahydropterin synthase/dimerization cofactor of hepatocyte nuclear factor 1 alpha (TCF1) 2; DcoH-like protein DCoHm; HNF-1-alpha		

For research use only.

GENE DESCRIPTION	pterin-4 alpha-carbinolamine dehydratase/dimerization cofactor of hepatocyte nuclear factor 1 alpha (TCF1) 2
GENE ID	84105
ORIGIN SPECIES	Human
ACCESSIONS	BC054021
VECTOR	pDONR223
NCBI SYMBOL	PCBD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggtgctcggggcgctcggggcgacgcggcgctgttggcggcgctcgaggccagagcctagggctagcggccatgtcatcagggtactcacagggtgactgcagaggagaggaaccaagctatactgacctaaagcagcaggatggtcggaattaagtgagagaga
tgccatctacaagaattctcctccacaattttaacaggcatttggctttatgtcccagttgcccacaagcagagaagatgaatcatcaccagaatgttcaatgtatacaacaagggtccagataactctcacctcacatgactgtggtgaactgacccaaaaagatgtgaagctggccaa
gtttattgaaaaagcagctgcttctgtgtga