

PAGE2-P antigen family, member 2 (prostate associated) Gene

Reference: PTXBC054022

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT16.4; GAGEC2; GAGEE2; PAGE-2; P antigen family member 2; G antigen family C 2; G antigen, family E, 2; P antigen family, member 2 (prostate associated); g antigen family E member 2; prostate-associated gene 2 protein; PAGE family member 2		

For research use only.

GENE DESCRIPTION	P antigen family, member 2 (prostate associated)
GENE ID	203569
ORIGIN SPECIES	Human
ACCESSIONS	BC054022
VECTOR	pDONR223
NCBI SYMBOL	PAGE2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgagcttctaagagcaagatccaatctcagaaagaggaaatgaccaagagcttccagccggttgatctgtgattgtccaggagccactgagggaaacgtcaagaagagggaaccaccaactgataatcagggtattgcacctagtggggagattgaaatcaagcagtg
ctgcttttcaaggccctgacatggaagcttttcaacaggaactggctctgcttaagatagaggatgagcctggagatggctctgatgicaggaggaggattatgccacttttgatctcactaaagtctggaagcagggtgatgcgaaccatag