

MRPS17-mitochondrial ribosomal protein S17 Gene

Reference: PTXBC054031

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC011; MRP-S17; RPMS17; S17mt; 28S ribosomal protein S17, mitochondrial; mitochondrial ribosomal protein S17		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S17
GENE ID	51373
ORIGIN SPECIES	Human
ACCESSIONS	BC054031
VECTOR	pDONR223
NCBI SYMBOL	MRPS17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgtagttcgtcatccgtccatgccagatggattgtgggaaggtgattgggacaaaatgcaaaagactgctaaagtgagagtgaccagctgttcttgatccctatttattaaagtattttaataagcggaaaacactcttgcacagatgcccttcagcagtgccacagttggggat
attgtgcttcagagctttacgtgtccacgagcaagcatgtgaacatgaactggtgagatcgttttcaaagttggaaaatcatagatccagtgacaggaaagccctgtgctggaactacctacctggagagtcggttgagttcggaaaccaccagctaagcaaaaatctggaagaa
ctcaatatctcttcagcacagtga