

TMEM69-transmembrane protein 69 Gene

Reference: PTXBC054041

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf154; transmembrane protein 69		

For research use only.

GENE DESCRIPTION	transmembrane protein 69
GENE ID	51249
ORIGIN SPECIES	Human
ACCESSIONS	BC054041
VECTOR	pDONR223
NCBI SYMBOL	TMEM69
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccaaggcccttgcttctctcatcattccagcgtatatgagcaagacacagtgctatcatatccccctgcagctttaa