

## RNF115-ring finger protein 115 Gene

Reference: PTXBC054049

|                |                                                                                                               |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                           | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                          | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | E3 ubiquitin-protein ligase RNF115; BCA2; ZNF364; rabring 7; zinc finger protein 364; ring finger protein 115 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | ring finger protein 115                                                                           |
| <b>GENE ID</b>          | 27246                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC054049                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | RNF115                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcggaggcttggcggccggggcgactcgggcgcgcgtgtagccgccaccgggttttctgccattttgcaaggcgaggctcagcccaactacgggaatatatgtccagatgtgaatcaggctttattgaagaagtacagatgattccagtttttagtggtggcggcagtcggatagacaataccacaacaacacattttgcagagctttggggccatttggatcacacgatgtttttcaagatttttagaccctttctaaagtagcagtcacactggaccaagataatagagccaatgaaaggggtcaccagactcacctgacttctggggagcaagacctccacgggtgccattgggtcggagatacagatctcaggaaagtctcgtctgcagatctccagctattgaaggaaactacaacacatctttgcaggatctttgcaaatctgcgaattctggatctccacaccttttctggagcgggatgctgcactccaacctggggactatgctgggtcagacagggcttgatgccattgtaaccagcttttaggacaactggaaaacacaggccctccccagctgacaaggaaaagatcacatctctccaacagtgcagtaactcagggaacaagttgatatgggttttagagtgccagatgcaagaagattacacagttgaagaggaaagtcggcagttaccttgcaatcacttttcacagcagttgattgtgcccgtggctagaactgcatgacatgtcctgtatgtaggaagagcttaaatggtgaggactctactcgcaaggccagagcactgaggcctctgcaagcaacagattagcaatgacagtcagctacatgaccgatggactttctga