

LYZL6-lysozyme-like 6 Gene

Reference: PTXBC054481

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-S-6a; LYC1; LYZB; PRO1485; TKAL754; UNQ754; lysozyme-like protein 6; epididymis secretory sperm binding protein Li 6a; lysozyme B; lysozyme homolog; lysozyme like 6		

For research use only.

GENE DESCRIPTION	lysozyme-like 6
GENE ID	57151
ORIGIN SPECIES	Human
ACCESSIONS	BC054481
VECTOR	pDONR223
NCBI SYMBOL	LYZL6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaaaggcgtactcatctatttgggtcagcagcttcttgcctaaatcagccagcctcatcgtcgtgtgacttggccaggtgctgcagctggaggacttgatgggtttgagggttactccctgagtgactggctgtgctgcttgggaaagcaagttcaacatacaagata
aafgaaatgcagacggaagctttgactatgacctttccagatcaacagccactactggtgcaacgattataagagttactcggaaaacctttgccacgtagactgtcaagatctgctgaatcccaaccttctgcaggcatccactgcgcaaaaaggattgtgtccggagcacgggggatg
aacaactgggtagaatggaggttgactgttcaggccggccactcttctactggctgacaggatgccgcctgagatga