

**PLEKHG4-pleckstrin homology domain containing, family G (with RhoGef domain)
member 4 Gene**

Reference: PTXBC054486

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARHGEF44; PRTPHN1; SCA4; puratrophin-1; PH domain-containing family G member 4; Purkinje cell atrophy associated protein 1; pleckstrin homology domain containing, family G (with RhoGef domain) member 4; pleckstrin homology and RhoGEF domain containing G4		

For research use only.

GENE DESCRIPTION	pleckstrin homology domain containing, family G (with RhoGef domain) member 4
GENE ID	25894
ORIGIN SPECIES	Human
ACCESSIONS	BC054486
VECTOR	pDONR223
NCBI SYMBOL	PLEKHG4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr