

MRPS24-mitochondrial ribosomal protein S24 Gene

Reference: PTXBC054865

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC335; MRP-S24; S24mt; bMRP-47; bMRP47; 28S ribosomal protein S24, mitochondrial; mitochondrial 28S ribosomal protein S24; mitochondrial ribosomal protein S24		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S24
GENE ID	64951
ORIGIN SPECIES	Human
ACCESSIONS	BC054865
VECTOR	pDONR223
NCBI SYMBOL	MRPS24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccctcgtgtgcagcgggttgcggggccacgggtgctgtcctggagccgagagctgcttgcgcttggcgccctgcacacctccccgggtcgcgcaagaacggcgccagagtacgcgtaagcaagggggacacgggtgacctacgaggaggcacacgcg
ccgcactacatcgccacccgtaaaaggctggctgtcgtgcacacaggtaacctggatggagaggacatgccgagagcgcaacagtgaggatgttttccttcgcaagttcatgtgggttaccttccaggctgcctggctgaccagctggttttaagcgccgggtaaccagttggag
atctgtgccgtgctcctgaggcagttgtctccacacaagtaacttctcgtgggtacagtgaaaccttgcctctacttttacaagtctcgtgcgactccacctccaactgtgccctcaaggttgtataagtaacctctag