

FXYP4-FXPYD domain containing ion transport regulator 4 Gene

Reference: PTXBC054876

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CHIF; FXYP domain-containing ion transport regulator 4; channel-inducing factor; FXYP domain containing ion transport regulator 4		

For research use only.

GENE DESCRIPTION	FXYP domain containing ion transport regulator 4
GENE ID	53828
ORIGIN SPECIES	Human
ACCESSIONS	BC054876
VECTOR	pDONR223
NCBI SYMBOL	FXYP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagagagtgcacctggcccttctctactggcagcctgactgccttggagccaatgacccattggccaataagacgaccccttctactatgactggaaaaacctgcagctgagcggactgatctgcggagggtcctggccattgctgggatcgccgagttctgagtgccaaatg
caaatgcaagagcagccagaagcagcacagtctgtacctgagaaggccatcccatcatcactccaggtctgcccactacttgcctga