

ANG-angiogenin, ribonuclease, RNase A family, 5 Gene

Reference: PTXBC054880

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALS9; HEL168; RAA1; RNASE4; RNASE5; angiogenin; RNase 5; angiogenin, ribonuclease, RNase A family, 5; epididymis luminal protein 168; ribonuclease 5; ribonuclease A A1; ribonuclease A family member 5		

For research use only.

GENE DESCRIPTION	angiogenin, ribonuclease, RNase A family, 5
GENE ID	283
ORIGIN SPECIES	Human
ACCESSIONS	BC054880
VECTOR	pDONR223
NCBI SYMBOL	ANG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgatggcctggcggtttgtgtgtgtctctgtctgggtctgggtctgacccaccgacccctggctcaggataactccaggtacacacacttctgaccagcactatgatgccaaccacagggccgggatgacagatactgtgaaagcatcatgaggagacggggcctgacct
cacctgcaaagacatcaacacatttattcatggcaaacgagcagcatcaggccatctgtgaaaacaagaatggaaccctcacagagaaaacctaagaataagcaagctcttcttcagggtcaccacttgcagctacatggaggttccccctggcctccatgccagtaccgagccac
agcggggttcagaaacgtgtgtgtgtgtgaaatggcttacctgtccacttgatcagtcgaatttccgtcgtccgtaa