

DYNLRB2-dynein, light chain, roadblock-type 2 Gene

Reference: PTXBC054892

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DNCL2B; DNLC2B; ROBLD2; dynein light chain roadblock-type 2; bithoraxoid-like protein; dynein light chain 2B, cytoplasmic; dynein, cytoplasmic, light polypeptide 2B; roadblock domain-containing protein 2; testicular tissue protein Li 53		

For research use only.

GENE DESCRIPTION	dynein, light chain, roadblock-type 2
GENE ID	83657
ORIGIN SPECIES	Human
ACCESSIONS	BC054892
VECTOR	pDONR223
NCBI SYMBOL	DYNLRB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaggtggaggaaaccttaagaggatccagatcataaaggggtattggaactatgggtgtaaatgcagaagatccacgtttgtgacctcaagtatgcaagacagctgactgccttcagagaaaaagtggccaacaaatcacaataaaagattag