

IGLV2-14-immunoglobulin lambda variable 2-14 Gene

Reference: PTXBC054893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	immunoglobulin lambda variable 2-14
GENE ID	28815
ORIGIN SPECIES	Human
ACCESSIONS	BC054893
VECTOR	pDONR223
NCBI SYMBOL	IGLV2-14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgtcttcaggccccctctcgcaggctgtgctgactcagccgtctcctcctcggcatctcctggagcatcagcagctcactgcaccctgcgcagaggctttatgtttatgactacaggatattgtgtaccaacagaagtcaggagggtctccccagtatctctgaggcacag
atcagatcagataaacaacagggtctctggagtcctccagccgtctctcggatccaagatgcttcggccaatgcaggaatttttagtcatctctgggtcctcggctcgaagatgaggctgactattactgtatgtttggcacaacagtccttgggtgttcggcggaggaccaggtgaccgtc
ctaagtacagcccaaggctgccccctcggctactctgttcccgccctcctcgtgaggtcctcaagccaacaaggccacactggtgtgtctcataaagtgacttctacccgggagccgtgacagtggcctggaaggcagatagcagccccgtcaaggcgggagtgagaccacacacccctc
caacaagaacacaagtacggcgccagcagctacctgagcctgacgctgagcagtggaagtcccaaaaagctacagctgccaggtcacgcatgaaggagcacctgtggagaagacagtgggccctacagaatgttcatag