

MRPL21-mitochondrial ribosomal protein L21 Gene

Reference: PTXBC055088

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L21mt; MRP-L21; 39S ribosomal protein L21, mitochondrial; mitochondrial ribosomal protein L21		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L21
GENE ID	219927
ORIGIN SPECIES	Human
ACCESSIONS	BC055088
VECTOR	pDONR223
NCBI SYMBOL	MRPL21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgccatggcagcatctccctgacggcgtacgttagggcggtggtcggtgcagccacagcatctgagacctcggggccggagcagcctcctttggtctgcttcgaaggtcaattcacagagcacttcatactaccaggatattgtctaaaacatccctgagt
tcaccaccttggccagaagttgttctgccagaccagttgaggagaccagacacacatgcagaggtcgtgaaagaaggtgaatgagatgatcgtcacggggcagtatggcaggctctttgccgtgggtgacatttgcagccgccagtggaagtgaccttgaagacctgacttaattggaa
atgaactagaccttgcgtgtggagagagaattcgactggagaaggtcctgctggttggggcagacaactcagcctgcttggcaagccactcctcgaaaggatctgttcgagtagaagccacagtcattgaaaagacagaatcatggccaagaatcattatgagattcagaaaaggaa
aaactcaagaagaaaagaatcgtcacgacccgcagactgtctccggataaacagcattgagattgctcgtgtttgttga