

EMG1-EMG1 nucleolar protein homolog (*S. cerevisiae*) Gene

Reference: PTXBC055314

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EMG1, N1-specific pseudouridine methyltransferase; EMG1 nucleolar protein homolog; C2F; Grcc2f; NEP1; ribosomal RNA small subunit methyltransferase NEP1; 18S rRNA (pseudouridine(1248)-N1)-methyltransferase; 18S rRNA (pseudouridine-N1)-methyltransferase N		

For research use only.

GENE DESCRIPTION	EMG1 nucleolar protein homolog (<i>S. cerevisiae</i>)
GENE ID	10436
ORIGIN SPECIES	Human
ACCESSIONS	BC055314
VECTOR	pDONR223
NCBI SYMBOL	EMG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcccagtgatggattcaagcctcgtgaacgaagcgggtggggagcaggcagcagctgggatgctctgccaccaagcgccccgactaggggcaggaaacaagatcgaggccgtaggcttattgtgtgctggaagggccagctggagacagtcaggtaggga
agacatatgagctactcaactgtgacaagcacaagctctatattgttgaagaatggacgggacctgggggaagcgcggccagatatcaccaccagagtttctgatgctgatggatagctccctgaaccgagctggcttgctacaggtttatccatcacagaagaatgttctgattgaagt
gaatccccagaccgaattcccagaacctttgaccgcttttggcctcatgtttcaacttttacacaagctcagtggttcgagcagctgatggccccagaagcctttgaaggtaattaagaatccagatcagatcacttccagttggatgatgaaagtggcacttctttccatcccgtgttc
agtgatgtcgtgagctgggtcccagcagtgatcctattgttttggtaggggctttgccatggcaaggtcagtggtggatatacagagaagatggtgtccatcagtaactaccccttctgctgccctcacctgtgcaaaactaccacagccttgagggaagtatgggggtcattga