

RHOQ-ras homolog gene family, member Q Gene

Reference: PTXBC056154

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rho-related GTP-binding protein RhoQ; ARHQ; HEL-S-42; RASL7A; TC10A; RAS-like, family 7, member A; epididymis secretory protein Li 42; ras homolog gene family, member Q; ras-like protein TC10; ras-like protein family member 7A; ras homolog family member Q		

For research use only.

GENE DESCRIPTION	ras homolog gene family, member Q
GENE ID	23433
ORIGIN SPECIES	Human
ACCESSIONS	BC056154
VECTOR	pDONR223
NCBI SYMBOL	RHOQ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcacggcccgccgctgatgctcaagtgcgtgggtgctcggcgacggggcggggcaagacgtgcctactcatgagctatgccaacgacgcctcccggaggagtagctgccaccgtcttcgaccactacgcagtcagcgtcaccgtgggggcaagcagtagctcctc
ggactctatgacacggccggacaggaagactatgaccgtctgagcccttacttacccaatgaccgatcttccttatatgctctcgtgtaaaatccagcctcattcaaaaatgtgaaagaggagtggtaccggaactaaggaaatagcaccaaatgtacccttttattaataggaactc
agattgatctccgagatgacccaaaacttagcaagactgaatgatataaagaaaaacatatgtgtggaacaaggacagaactagcaaaagagataggagcatgctgctatgggaatgttcagcttaaccagaagggaatgaagactgtttgatgaggctatagccattta
actccaaagaacacactgtaaaaaaagaatagatcaagatgataaactgttttaattacgtga