



Reference: PTXBC056219

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP-dependent RNA helicase DHX29; DDX29; DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide 29; DEAH (Asp-Glu-Ala-His) box polypeptide 29; DEAH box protein 29; DEAH-box helicase 29; nucleic acid helicase DDXx; DExH-box helicase 29		

For research use only.

GENE DESCRIPTION	DEAH (Asp-Glu-Ala-His) box polypeptide 29
GENE ID	54505
ORIGIN SPECIES	Human
ACCESSIONS	BC056219
VECTOR	pDONR223
NCBI SYMBOL	DHX29
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]