

SPRR1B-small proline-rich protein 1B (cornifin) Gene

Reference: PTXBC056240

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CORNIFIN; GADD33; SPR-IB; SPRR1; cornifin-B; 14.9 kDa pancornulin; small proline-rich protein 1B; small proline rich protein 1B		

For research use only.

GENE DESCRIPTION	small proline-rich protein 1B (cornifin)
GENE ID	6699
ORIGIN SPECIES	Human
ACCESSIONS	BC056240
VECTOR	pDONR223
NCBI SYMBOL	SPRR1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttccagcagcagagcagccctgcatccacccctcagcttcagcagcagcaggtgaaacagcctgccagcctccacctcaggaaccatgcaccccaaaacaaaggagccctgccaccccaaggtgcctgagccctgccaccccaagtcctgagccctgccagccc
aagggtccagagccatgccaccccaaggtgcctgagccctgccctcaatagtcactccagcaccagccagcagagacccaagcagaagtaa