

ZNF254-zinc finger protein 254 Gene

Reference: PTXBC056243

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BMZF-5; HD-ZNF1; ZNF539; ZNF91L; zinc finger protein 254; CTD-2017D11.1; bone marrow zinc finger 5; hematopoietic cell-derived zinc finger protein 1; zinc finger protein 539		

For research use only.

GENE DESCRIPTION	zinc finger protein 254
GENE ID	9534
ORIGIN SPECIES	Human
ACCESSIONS	BC056243
VECTOR	pDONR223
NCBI SYMBOL	ZNF254
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaggacccctagaaagcctagaaatgggactgttgacatttagggatgtggccatagaattctctctggaggagtggaacacctggacattgcacagcagaatttatagaaatgtgatgttagagaactacagaacacctggccttctgggtattgctgtctctaaagccagacctgat
cacctgtctggaacaagggaagagccctggaatatgaagcgacatgagatggatgaaccccgaggattgatttttactactgtga