

LEPROT-leptin receptor overlapping transcript Gene

Reference: PTXBC056250

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LEPR; OB-RGRP; OBRGRP; VPS55; leptin receptor gene-related protein; DnaJ (Hsp40) homolog, subfamily C, member 6; OB-R gene-related protein; endospanin-1; leptin receptor overlapping transcript protein; leptin receptor overlapping transcript		

For research use only.

GENE DESCRIPTION	leptin receptor overlapping transcript
GENE ID	54741
ORIGIN SPECIES	Human
ACCESSIONS	BC056250
VECTOR	pDONR223
NCBI SYMBOL	LEPROT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggttaaagctctcgtggcattatcctcagtgggcctattggactgactttcttatgctgggatgtgccttagaggattatggcgtttactggccttattcgtctgattttccacgccatctccccatccccattcattgccaagagtcacctgactcagatgcaaccagta
gtgcctgtcgggaactggcattatctcactactggaattggtttctgctttggattcctgtattctgtctgtgctgtgatcaaatggggagcctgcccctgtgtggcaggcaatgcatcattttcctacaatcaagggttttctctatatttgaagaggagatgatttagctgg
gagcagtggtag