

AP3M2-adaptor-related protein complex 3, mu 2 subunit Gene

Reference: PTXBC056257

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AP47B; CLA20; P47B; AP-3 complex subunit mu-2; HA1 47 kDa subunit homolog 2; HA1 47kDa subunit homolog 2; adapter-related protein complex 3 mu-2 subunit; adapter-related protein complex 3 subunit mu-2; adaptor-related protein complex 3 subunit mu-2; clath		

For research use only.

GENE DESCRIPTION	adaptor-related protein complex 3, mu 2 subunit
GENE ID	10947
ORIGIN SPECIES	Human
ACCESSIONS	BC056257
VECTOR	pDONR223
NCBI SYMBOL	AP3M2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatccatagtcctttgatcaactcctctggagacatttctggagaaacattggaaaagtgtggtcagccgttctgttctgattactttttgagggcgaagagagagctactaggcagaaaatgtgcctccggttatccctaccctcaccactatctttaagtgtttaccgccacaaga
tctttttgtggccgtgatccagacggagggtcccccctctgtttgtcattgagtttcttaccgagtggtggacacatttcaggattattttgagtgctgttcagagccagtgatcaagacaatgtagttgtgtttatgaggtaattggaagagatgcttgacaatggtttccattggctaccgagtcg
aacattcttaagaactcataaagcctcctaccatccttgaacggtttgcaacacatcacaggaagcagcaatgtgggtgaccagctccactggcgagctgtcagtggtgcttggcgacgagctgggtgaaatataccaacaatgagggcctattttgatgtgattgaagagattgatg
caattattgataaatcaggtccacaattactgctgagatccaggggtgattgacctgtgtcaagctgactgcatgccagacctacatttcttcataaacctaggtgttgatgatgcagcttccatcttgtgttcgtttcaaacgctgggaatctgagcgcacccctctctcatccc
tctgatggaaactccgctgctgttcttaccatgtcagtgacagaaatgctgtcttgggatgtag