

## C11orf84-chromosome 11 open reading frame 84 Gene

Reference: PTXBC056402

|                               |                                                                       |                     |                |
|-------------------------------|-----------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                   | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                  | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | uncharacterized protein C11orf84; chromosome 11 open reading frame 84 |                     |                |
| <b>For research use only.</b> |                                                                       |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 11 open reading frame 84                                                               |
| <b>GENE ID</b>          | 144097                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC056402                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C11orf84                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggccctaaaggccgagggcgccgactcactgcttcgaggtgacgtgaaatgcgaggaaggggagggagggagggccatggtggcgccgtaattccggcccgagccgatgctcagagtgacccaacaggagaagacccacccgctagaccagcccctaga  
ggcagggcagtgatggctgtgaggagccgaagcagcaggtgtcttgggagcaggagttcctgtgggagcagcccagggagggcggcggcactgtgcatgtgtgtgtggcgtgagatccgggacccctcggccgacacagctcgtcgcacatcttgagcagcaccctcac  
acctggacgtgagccctctgagaagagcaatctcggagggcctggagtgaaaggggtggccctcttgcagacgtgagagctgagcagccgtcccccacccaactcagactcgggcccaggtgcccacccagaccagacccaacccagacgtgccagaatgccagccgaat  
cgtggttctccttgacttgaggataacccatccctccctaaagagcagggcccaggggactccgcccctcagcttctcgtctccctgccacagagccaggaataagaagccccgtgtgagagatggaaggaaacccagggggaagagccagtcagaaagaaaagagga  
gacatgacacaaaacgtgacccctgaccagagccccatgccagactcggccacggagacttgcagcaccagccgaggtccgacactcactgacggcagcttccccggtctcgttgcagctcttccccacacccagctcagggggccagacagcaaggactcacc  
caagacagggaggtggcagaagggagccttccccggcgagagccccctcagctccccctccgggctccgagacactggatctccaggttatccgctgctggatggaggagccccagcggctcagctcctcgaagactgtccaggcaccacccagggcaccagcg  
tgtgggagcaggtgacacctcagactggccacagttgtgcagaatccagcaccactgtggcagggaagcgggaaaaagggaatggagtgtaa