

**NSUN5C-NOL1/NOP2/Sun domain family, member 5C Gene**

Reference: PTXBC056405

|                |      |                     |                |
|----------------|------|---------------------|----------------|
| Product type:  | DNA  | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       |      |                     |                |

**For research use only.**

---

|                         |   |
|-------------------------|---|
| <b>GENE DESCRIPTION</b> | NOL1/NOP2/Sun domain family, member 5C  |
| <b>GENE ID</b>          | 260294  |
| <b>ORIGIN SPECIES</b>   | Human   |
| <b>ACCESSIONS</b>       | BC056405  |
| <b>VECTOR</b>           | pDONR223  |
| <b>NCBI SYMBOL</b>      | NSUN5C  |
| <b>LENGTH</b>           | ProteoGenix bp  |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

atgccgagcagacagctggaggagccggggcaggggacacctagcccgtgctgctgcatgccctggcagggtccagcagcagccctgtccacggctcacttccctccctgcagcggctctactccatgtgctccctctgccaggagagaatgaagacatggtacaag  
atgcgctgcagcagaacccggcgccctcagctagctcccgcctgctgcccggccccaccaggcctgagcagcttcccgggtgccgagcactgcctccggctcccccaagaccacgcttagcgggtgcttctctgtgctgtaattgaacgggtcgagatgccgacgtga