

CROP-cisplatin resistance-associated overexpressed protein Gene

Reference: PTXBC056409

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CROP; CRA; CREAP-1; LUC7A; OA48-18; hLuc7A; luc7-like protein 3; CRE-associated protein 1; LUC7-like 3; cAMP regulatory element-associated protein 1; cisplatin resistance-associated-overexpressed protein; okadaic acid-inducible phosphoprotein OA48-18; LUC		

For research use only.

GENE DESCRIPTION	cisplatin resistance-associated overexpressed protein
GENE ID	51747
ORIGIN SPECIES	Human
ACCESSIONS	BC056409
VECTOR	pDONR223
NCBI SYMBOL	CROP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatttcgcccgcagttgttgatgagttatggccgggaccgaaacctagccccggacgagaagcgcagcaacgtgcggtgggaccacgagagcgtttgtaaatattatctctgtgtttttgtcctgcggaattgtcacaatacacgttctgatcttgatgtattggaaggagaga
taacattagtgatgtcagcaaatgtttggaagatgacaagtggatggaggagtag